


Buena calidad ajustable espalda corrector de postura corrector de soporte de postura, banda de corrección jorobado simple

Nombre	Corrector de postura / Corrector de postura
Uso	Prevenga el encorvamiento, mejore su postura y alivie el dolor Su espalda y hombros se ajustarán lentamente y desarrollará memoria muscular para la postura firme
tamaño	Talla única
Tamaño de cofre adecuado	75-115
Logotipo Personalizar	Personalizar el logotipo de los clientes

## SUITABLE CHEST SIZE



Pluma

### FOAM PERFORATED FABRIC

soft and thin, moderate thickness,  
full material ouch design.



### MIDDLE SUPPORT OF

Using positive pu-skin, back  
with sponge, comfortable to  
wear.



### FIXED SIDE BUCKLE SUPPORT

Square buckle position just designed at both  
ends of the shoulder strap, will not cause a  
bit uncomfortable.



# How to use

1.



Place posture corrector brace, it's in a "Figure 8 shaped"

2.



Slip the brace on like you would put on a backpack. The triangle buckle should point down once the brace is on your back.

3.



To adjust the fit, undo the velcro closures and pull the straps tighter until your desired level of posture support is reached.



You should feel your shoulders being gently pulled back.